


Equip™

Standard Blade Range

*Gold Standard for Value
Gold Standard for Quality*



Equip™ standard blades bring new meaning to the "Gold Standard" in laryngoscopes.

Value for money without compromise on quality.

New streamlined profile.

Conforms to ISO 7376-1 and ASTM F 965-99
Approvals - ISO 9001, ISO 13485
FDA registered manufacturer



 **TRUPHATEK®**

Technical Specifications



Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.	Krypton Bulb
Mac0	Neo-natal	77	9	5810/	Small
Mac1	Infants	90	9	5811/	Small
Mac2	Infants	100	12	5812/	Small
Mac3	Adults	130	15	5813/	Large
Mac4	Adults	155	15	5814/	Large
Mac5	Adults	172	15	5815/	Large



Size	Finish	Batteries	Order No.
UltraSafe™ Standard	Knurled	C	5525
UltraSafe™ Mini	Knurled	AA	5526
UltraSafe™ Paediatric	Knurled	N	5527
Stubby	Knurled	AA	5523
Dolphin™ Standard	Knurled	C	5571
Dolphin™ Mini	Knurled	AA	5572



Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.	Krypton Bulb
Miller00	Neo-natal	64	11.5	5829/	Small
Miller0	Neo-natal	75	11.5	5830/	Small
Miller1	Infants	102	11.5	5831/	Small
Miller2	Infants	155	14	5832/	Large
Miller3	Adults	195	14	5833/	Large
Miller4	Adults	205	14	5834/	Large



Bulb Specifications					
Description	Screw	Amp.	Volt	Diameter	Order No.
Krypton Small	1/8 - 72 NS	0.3A	2.5V	4.0 mm	512201
Krypton Large	8 - 32 UNC	0.3A	2.5V	5.4 mm	512202



Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.	Krypton Bulb
Oxy00	Neo-natal	76	11.5	5827/	Small
Oxy0	Neo-natal	87	11.5	5828/	Small



Boxed sets	Order No.
Handle Plus 3 Blades Mac 2,3,4	5713E/
Handle Plus 4 Blades Mac 1,2,3,4	5711E/
Handle Plus 3 Blades Miller 2,3,4	5717E/

* Additional combinations available upon request

Truphatek®, Dolphin™, UltraSafe™ and Equip™ are trademarks of Truphatek International and should not be used without prior written permission from Truphatek.

From Truphatek's Warranty Statement:

"All Truphatek products are unconditionally warranted to be of good quality and free from manufacturing defects in material and workmanship for a period of 5 years from date of manufacture. Replacement product will be provided without charge for any Truphatek product not meeting this warranty."

Full text of Truphatek's warranty can be viewed at our website www.truphatek.com or by request from Customer Service Desk at helpdsk@truphatek.com

Distributed by

TRUPHATEK
www.truphatek.com

Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel
Tel: 972-9-8851155, Fax: 972-9-8851212, helpdsk@truphatek.com

© 2005-2007 TRUPHATEK

EQ-1-02/07 Rev. 3