

GreenLite™

Disposable Metal Laryngoscope Blade

*Single Use Disposable
Made of Steel!*



*Uncompromising Performance in
the fight against cross infection*

- *Unique patented design duplicates metal blade performance and feel*
- *Excellent light transmission*
- *Provides visualization without "dazzle" effect*
- *Simple, reliable "click" locking and unlocking mechanism*
- *Fits onto Green system handles*

Conforms to ISO 7376-3 and ASTM F 1195-99
Approvals - ISO 9001, ISO 13485
FDA registered manufacturer
Patent Pending



 **TRUPHATEK®**

Technical Specifications



Order No.	Description	Overall Length (mm)	Distal Width (mm)
4341	Mac 1	92	9
4342	Mac 2	102	12
4343	Mac 3	132	15
4344	Mac 4	157	15



Order No.	Description	Overall Length (mm)	Distal Width (mm)
4319	Miller 00	68	11.5
4320	Miller 0	78	11.5
4321	Miller 1	105	11.5
4322	Miller 2	158	14
4323	Miller 3	198	14



Order No.	Description	Comments
4306	EMS pack	Handle Plus 3 blades - Mac 2,3 & 4

From Truphatek's Warranty Statement:

"GreenLite™ in original packaging is unconditionally guaranteed to be of good quality and free from manufacturing defects in material and workmanship for a period of 3 years from date of manufacture. Replacement product will be provided without charge for any Truphatek product not meeting this warranty."

Full text of Truphatek's warranty can be viewed at our website www.truphatek.com or by request from Customer Services at helpdsk@truphatek.com

Truphatek™ and GreenLite™ are trademarks of Truphatek International Ltd. and are not to be used without prior written permission of the manufacturer.

Distributed by



Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel
Tel: 972-9-8851155, Fax: 972-9-8851212, helpdsk@truphatek.com