

**Green**  
spec-2™

*The Most Extensive Range of  
Fibre Optic Laryngoscope Blades  
with the new FibreClip™  
Replaceable Light Guide*



- *Widest choice of Macintosh and Miller Blades and the unique Oxyport paediatric blades*
- *Incomparable cold, bright light - 2800 LUX*
- *Simple, trouble-free use. No screws, no lost parts*
- *Autoclavable FibreClip™ light guide, tested to withstand over 500 cycles*
- *Stainless steel parts guaranteed for 5 years*

Approvals - ISO 9001, ISO 13485  
Conforms to ISO 7376-3 and  
ASTM F 1195-99  
FDA registered manufacturer



**TRUPHATEK®**

# Technical Specifications

Blades				
Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.
Mac0	Neo-natal	77	9	4230/
Mac1	Neo-natal	90	9	4231/
Mac2	Infant	100	12	4232/
Mac3	Adult	130	15	4233/
Mac4	Adult	155	15	4234/
Mac5	Adult	172	15	4235/

Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.
Miller00	Neo-natal	64	11.5	4249/
Miller0	Neo-natal	75	11.5	4250/
Miller1	Infant	102	11.5	4251/
Miller2	Child	155	14	4252/
Miller3	Adult	195	14	4253/
Miller4	Adult	205	14	4254/

Size	Use	Overall Length (mm)	Distal Width (mm)	Order No.
OxyMiller00	Neo-natal	76	11.5	4240/
OxyMiller0	Neo-natal	87	11.5	4241/

Sets Include*		
Green SpeX™ Handle	GreenSpec 2™ Blades with FibreClip™ bundle	Order No.
Standard size	Mac 1,2 & 3	4260/X
Standard size	Mac 1,2,3 & 4	4261/X
Standard size	Mac 2,3 & 4	4262/X
Standard size	Miller 00,1, 2 & 3	4263/X

\* Additional combinations available upon request

TruCharge™ Power Station	Order No.
Charger base, 2 Green SpeX-R™ Rechargeable handles and 2 Rechargeable Lithium-Ion Battery packs	5035 E, U/J, B
Charger base, 1 Green SpeX-R™ Medium Rechargeable handle, 1 Green SpeX-R™ Mini Rechargeable handle and 2 Rechargeable Lithium-Ion Battery packs	5036 E, U/J, B

Bulbs	Order No. Box of 6
Green SpeX™ Handle Xenon 2.5V 0.7A 4.8 O.D. M3.5x0.35	512501/6
Green SpeX-R™ Rechargeable Handle Xenon 3.5V 0.7A 4.8 O.D. M3.5x0.35	513501/6

FibreClip™ Light Guide			
Size	Use	Overall Length (mm)	Order No. Pack of 2 FibreClip™ Light Guides
Mac0	Neo-natal	55	423002
Mac1	Neo-natal	55	423102
Mac2	Infant	67	423202
Mac3	Adult	91	423302
Mac4	Adult	117	423402
Mac5	Adult	117	423502
Miller00	Neo-natal	37	424902
Miller0	Neo-natal	43	425002
Miller1	Infant	71	425102
Miller2	Child	110	425202
Miller3	Adult	147	425302
Miller4	Adult	157	425402

GreenSpeX™ Handles			
Size	Finish	Batteries	Code
Standard	Knurled	C	4211X
Mini	Knurled	AA	4212X
Stubby	Knurled	AA	4213X

GreenSpeX-R™ Rechargeable Handles		
Size		Code
Standard	Medium Battery pack	4211X-R
Mini	Mini Battery pack	4212X-R

Rechargeable Battery Packs	
Size	Code
Medium	5071
Mini	5072

## From Truphatek®'s Warranty Statement:

"All Truphatek® products are unconditionally guaranteed to be of good quality and free from manufacturing defects in material and workmanship for a period of 5 years from date of manufacture. Replacement product will be provided without charge for any Truphatek® product not meeting this warranty. Fibre Clip™ replaceable fibre-optic light guides are guaranteed to perform for a period of 15 months from date of manufacture"

Full text of Truphatek®'s warranty can be viewed at our website [www.truphatek.com](http://www.truphatek.com) or by request from Customer Service Desk at [helpdesk@truphatek.com](mailto:helpdesk@truphatek.com)

Green Spec2® is a patented product. Green Spec®, Green SpeX™, Green SpeX-R™ and Green Spec2® are trademarks of Truphatek® and should not be used without prior written permission from Truphatek International. Ltd.

Distributed by

**TRUPHATEK**  
www.truphatek.com

Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel  
Tel: 972-9-8851155, Fax: 972-9-8851212, [helpdesk@truphatek.com](mailto:helpdesk@truphatek.com)

© 2005-2007 TRUPHATEK® All rights reserved

GS2-1-04/07 REV 05