



From Truphatek®, the new source for 'value-for-money' fibre optic cables, connectors and LED light sources.




- *Top Quality A2 fibre optic raw material from SCHOTT Germany providing high light transmission, large aperture angle and low colour shift*
- *High level polishing and stringent quality control*
- *Twelve month warranty with 'easy to read' manufacturing date*
- *First-class, locally available customer service*
- *OEM & Development Services available.*

FDA No. K063309
FDA registered manufacturer
Approvals - ISO 9001, ISO 13485



Technical Specifications



Instrument Connectors

Catalogue #	Fits	
65IN202AM	ACMI Male	
65IN206WDF	Wolf Female Snap On	
65IN207WDM	Wolf Male	
65IN201STO	Karl Storz	
65IN203AF	ACMI Snap On Female	
65IN205AC	ACMI / Codeman Clip-on	
65IN204SK	Stryker	

Light source connectors

Catalogue #	Fits	
65LS106ACD	ACMI / COD MAN	
65LS102ST	Storz	
65LS105WD	Extended Wolf	
65LS109FN	Fijinon	
65LS103PA	Pentax	
65LS101O	Olympus	
65LS104M	Machida	
65LS107F	Fuji	

Colour

Code	Options
GR	Grey
BL	Blue

Bundle Diameter

Code	Options
D3	3.5mm
D5	5.0mm

Length

Code	Options
L6	6 ft (1.8 m)
L7	7.5 ft (2.3 m)
L10	10 ft (3 m)

From Truphatek®'s Warranty Statement:
 Truphatek®'s Tru-Cable™ fibre optic cables and connectors carry a 12 month warranty from the date of manufacture against defects in manufacture and materials of the strain relief and connectors. Truphatek® does not warranty cable connectors against any damage caused by use with inappropriate light sources. Truphatek® does not warranty the cable against fibre breakage due to abuse or inappropriate handling. Full text of the Truphatek® warranty and further information about Truphatek®'s products are available on our website www.truphatek.com or through our Customer Service Desk at helpdesk@truphatek.com.

Note : Truphatek® and Tru-Cable™ are trademarks of Truphatek International Limited and are not to be used without prior written permission.

Distributed by

 **TRUPHATEK**
 www.truphatek.com

Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel
 Tel: 972-9-8851155, Fax: 972-9-8851212, helpdesk@truphatek.com

© 2007 TRUPHATEK All rights reserved

CAB-1-05/07