

## New Product Announcement – TruLite Intubation Kit

### *Full-Sized Intubation Kit for the TruLite Laryngoscope*



The new TruLite Intubation Kit is a full-sized deployment system for the innovative TruLite Disposable Laryngoscope. It holds 6 different TruLite Laryngoscopes, and everything else needed for ALS airway care.

The Intubation Kit has color-coded pockets for quickly differentiating Macs from Millers. It has clear window pockets for labels. It also features staggered pockets for ET tubes, and a clear large pocket for percutaneous device or lighted stylets. It also offers elastic loops for oral airways, magill forceps, and can even carry a backup specialty laryngoscope.



It is designed specifically for the lightest weight and brightest disposable steel-bladed laryngoscope in the world.

The [TruLite Laryngoscope](#) is a completely disposable, integrated handle and blade combination. The blade and handle are permanently attached together. Every blade has its own power supply from 2 AA batteries in its own compact and lightweight handle. They are offered in color-coded Mac 2, 3, and 4 sizes, and soon in Miller 1, 2, and 3 sizes.

- **Forged Stainless Steel Blades**
- **Super-Bright, blue-white LED illumination**
- **No Cleaning Cost**



The TruLite Intubation kit is manufactured from rugged waterproof; 500-denier black nylon polyurethane coated fabric and features a durable double-pull YKK zipper. Crafted to exacting standards by Propak Manufacturing, a world leader in premium quality cases and packs. Dimensions - H: 9.5" (24cm), L: 14" (35.6cm), D: 3" (7.62cm)

# 12621 Case Only Price \$69.95