

Pro™

Articulating Laryngoscope Blades

*“Proven clinically to provide greater flexibility and improved control in a wide range of difficult intubation cases.”**



- Incorporates a movable tip controlled by a lever on the handle of the Pro blade.
- Depressing the lever exposes the larynx and glottic opening and simplifies the elevation of the epiglottis.
- Unique anti-seize spring design ensures safety in use.
- Force applied during intubation is greatly reduced*
- Indicated in cases of c-spine precautions, morbid obesity and other causes of anterior glottic exposure*.

* J Accid Emerg Med 2000; 17:364-365
Can J of Anes. 48; 212-213 (2001)

Shucman Pro conforms to ISO 7376-1 and ASTM F 965-99
Green Spec Pro conforms ISO 7376-3 and ASTM F 1195-99
Approvals - ISO 9001, ISO 13485
FDA registered manufacturer



TRUPHATEK®

Technical Specifications

Shucman-Pro

Sets Contain:	Code for Set	Overall Blade Length (mm)	Distal Blade Width (mm)
Shucman Pro Blade size Mac 3.5 Dolphin™ Waterproof Handle 1 spare bulb	5770	138	15
Shucman Pro Blade size Mac 4.5 Dolphin™ Waterproof Handle 1 spare bulb	5771	159	15

Bulbs for Shucman Pro™ Blades	Code
Halogen - box of 6	512717/6

Green specpro

Sets Contain:	Code for Set	Overall Blade Length (mm)	Distal Blade Width (mm)
GreenSpec Pro Blade size Mac 3.5 GreenSpeX™ Handle 1 spare bulb	4770X	138	15
GreenSpec Pro Blade size Mac 4.5 GreenSpeX™ Handle 1 spare bulb	4774X	159	15

Bulbs for GreenSpeX™ Handles		Code
Xenon - box of 6	Xenon 2.5V 0.7A 4.8 O.D. M3.5x0.35	512501/6

From Truphatek's® Warranty Statement:

"All Truphatek® products are unconditionally guaranteed to be of good quality and free from manufacturing defects in material and workmanship for a period of 5 years from date of manufacture. Replacement product will be provided without charge for any Truphatek® product not meeting this warranty."

Full text of Truphatek's® warranty can be viewed at our website www.truphatek.com or by request from Customer Services at helpdsk@truphatek.com

Shucman Pro™, GreenSpec Pro™, Dolphin™ and GreenSpeX™ are trademarks of Truphatek International and should not be used without prior written permission from Truphatek International Ltd., Israel.

Distributed by

TRUPHATEK®
www.truphatek.com

Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel
Tel: 972-9-8851155, Fax: 972-9-8851212, helpdsk@truphatek.com

© 2005-2007 TRUPHATEK® All rights reserved PGSSH-01-09/07 Rev. 03