

TruCharge™

Power Station

Charger Base with Accessories



Economical and Environmentally Friendly Solution for charging laryngoscope handle batteries!

Custom designed charger base for use exclusively with Green SpeX-R™ laryngoscope handles and rechargeable Lithium Ion Battery Pack.

- Modern Lithium Ion technology:
Fastest charging
Longest battery life
No memory build up
- Safe, Reliable solution:
Multicolour LEDs show change of status
PCM element protects against short circuiting, overloading and over-discharging

Approvals - ISO 9001, ISO 13485
EN 60601, FDA registered manufacturer



TRUPHATEK®

Frequently Asked Questions

Why should my department change to rechargeable batteries?

- Cost Saving, Economical alternative
"...A single rechargeable... battery can replace up to 1000 single-use alkaline batteries over its lifetime"*
- Green Environment - With fewer batteries to dispose of you will be contributing to a cleaner and friendlier environment for all of us.

Why should my department prefer Lithium Ion rechargeable batteries?

Remember the performance advantages of Lithium Ion batteries:

- Lightweight
- Retain their charge, are not prone to the memory effect
- Work over a broad temperature range
- High voltage supports Xenon bulb with 50% higher light intensity

How is the rechargeable handle different?

How can I identify a rechargeable handle?

- Rechargeable handle is clearly marked as 'Green SpeX-R'
- Bottom of the handle has a black insert with a hole in the middle
- The bulb is Xenon 3.5 Volt

Why are there multicolored LEDs on the charger base?

- To ensure Reliability. Status of the handle is clearly visible. Red LED indicates charging is in process
Green LED indicates charging completed and assures you the battery will perform reliably
- To ensure Safety. Faulty, unreliable batteries will be clearly identified by the blinking red LED. A Battery Pack that has not completed the charging within 6 hours will trigger the faulty 'blinking red' LED sign and should be replaced

How much time does it take to charge a handle?

This can generally take anywhere from one to four hours. TruCharge™ Rechargeable Lithium Ion Battery Pack is supplied partially charged and should be fully charged prior to use.

How much time or how many intubations will one charge last for?

The capacity of the Battery Pack provides at least 1 hour of continuous work with Xenon bulb 3.5V 0.7A. The number of intubations depends on the length of each intubation, whether the Battery Pack was completely charged when work began and whether it is recharged between intubations

What is the life span of a Battery Pack?

TruCharge Rechargeable Lithium Ion Battery Pack is under warranty for one year from date of purchase or up to 500 charging cycles.

*Quotes from "INFORM: Strategies for a better environment"
Product no. 5071
Lot no. 07/05

Technical Specifications



Legend: E = EU plug, A = Australian plug, B = British plug, U/J = USA/Japanese plug

KITS	Ordering Code	Description
	5035 E, A, B or U/J	TruCharge™ Desk-Top Charger base with power supply and plug plus 2 medium size handles with 2 battery packs
	5036 E, A, B or U/J	TruCharge™ Desk-Top Charger base with power supply and plug plus 1 medium size and 1 mini size handle and 2 battery packs



BASES	Ordering Code	Description
	5020 E,A,B or U/J	TruCharge Desk-Top Charger base for 2 medium handles with power supply and plug
	5030 E, A, B or U/J	TruCharge Desk-Top Charger base for 1 medium and 1 mini handle with power supply and plug



BATTERIES	Ordering Code	Description
	5071	TruCharge™ Lithium-ion Rechargeable Battery Pack – medium size
	5072	TruCharge™ Lithium-ion Rechargeable Battery Pack – mini size



HANDLES	Ordering Code	Description
	4211XR	Green SpeX-R™ Rechargeable Laryngoscope Handle with Xenon Bulb – medium size. Battery not included
	4212XR	Green SpeX-R™ Rechargeable Laryngoscope Handle with Xenon Bulb – mini size. Battery not included



BULBS	Ordering Code	Description
	513501/6	Xenon Bulb 3.5V – Box of 6

From Truphatek's Warranty Statement:

We guarantee the TruCharge™ Charger Base and the TruCharge™ Rechargeable Battery Packs against faults or defects in material, design or workmanship for a period of one year from date of purchase provided the products are used in accordance with the instructions for use, only Truphatek® spare parts are used and no unauthorized modifications or repairs were made to the units.

Full text of Truphatek's warranty can be viewed at our website www.truphatek.com or by request from Customer Services at helpdsk@truphatek.com

Truphatek®, TruCharge™, Green SpeX-R™ are trademarks belonging to Truphatek International Ltd. and should not be used without the prior written permission of the owner.

TRUPHATEK
www.truphatek.com

Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel
Tel: 972-9-8851155, Fax: 972-9-8851212, helpdsk@truphatek.com

Distributed by