

# Lite-Blade™ Slims

Disposable Laryngoscope Blade

**Lite-Blade**  
is "...the best single-use  
plastic blade..."\*

Experts and Clinicians Agree



## Lite-Blade Slims™ Patented, Disposable Laryngoscope Blade

- Ergonomic slim design with tooth protection area
- Feels and performs exactly like a metal blade
- Disposable and "single-patient-use" blades decrease the high risk of cross-infection
- Tested for tensile strength well in excess of the physical forces required during intubation
- Unbeatable, strong, brilliant white light
- Individually packaged - ready for use
- Available in all Macintosh and Miller sizes

\*S. Rassam et al *Anaesthesia* 2005; 60: 384-394

\*S.J. Twigg et al *Br J Anaesth* 2003; 90: 8-13

Compatible with standard laryngoscope handles  
conforming to ISO 7376-1 & ASTM F 965-99  
IL patent no. 107,594 and associated patents registered worldwide  
Approvals - ISO 9001, ISO 13485  
FDA registered manufacturer




 **TRUPHATEK®**

# Technical Specifications



*Lite-Blade Slims*

Order No.	Description	Overall Length (mm)	Distal Width (mm)
4041	Mac 1	90	12
4042	Mac 2	106	13
4043	Mac 3	132	17.5
4044	Mac 4	156	17.5



*Lite-Blade*

Order No.	Description	Overall Length (mm)	Distal Width (mm)
4020	Miller 0	76	11.5
4021	Miller 1	103	11.5
4022	Miller 2	155	16
4023	Miller 3	195	16



Order No.	Description	Comments
4006	EMS pack	Handle Plus 3 blades - Mac 2,3 & 4
4005	EMS pack	Empty pouch

### From Truphatek®'s Warranty Statement:

"Lite-Blades™ in original packaging are unconditionally guaranteed to be of good quality and free from manufacturing defects in material and workmanship for a period of 5 years from date of manufacture. Replacement product will be provided without charge for product not meeting this warranty."

Full text of Truphatek®'s warranty can be viewed at our website [www.truphatek.com](http://www.truphatek.com) or by request from Customer Services at [helpdsk@truphatek.com](mailto:helpdsk@truphatek.com)

Lite-Blade™ is a patented product. Lite-Blades™ and Lite-Blade Slims™ are trademarks of Truphatek® and should not be used without prior written permission from Truphatek International. Ltd.

**Distributed by**



**Truphatek International Ltd. P.O.Box 8051 Netanya 42504, Israel**  
**Tel: 972-9-8851155, Fax: 972-9-8851212, [helpdsk@truphatek.com](mailto:helpdsk@truphatek.com)**

© 2002-2007 TRUPHATEK® All rights reserved

LBSL-1-05/07 Rev. 04